Thermoplastic material for denture bases

Product characteristics: Thermoplastic material in the form of a granulated product for making dentures using the pressure injection method, intended for use in dental and technical laboratories. The product does not contain methyl methacrylate.

Classification according to EN ISO 20795-1: Type 3**Chemical composition**

Granulated product: micromeritic polyamide, pigments

Performance:

Technical data necessary for correct material processing are given in Table 1.

Intended use of the product

Indications: total dentures (edentulous), partial, retention elements in framework dentures, micro dentures, splints and space maintainers.

Contraindications: Do not use the product in patients who are allergic to any of the ingredients.

Target group or group of target patients and intended users: People of different ages, with partial or complete lack of teeth, eligible for the prosthetic treatment with removable artificial teeth.

Side effects: Not found in the current period of the use of the product. In the absence of proper hygiene, irritations and oral inflammations caused by fungi and bacteria adjacent to the denture are possible.

Clinical benefits: Restoring the chewing function, correct phonetics and improving facial aesthetics, which positively affects the patient's psychological comfort.

Directions for use

Preparation of the working model: The working model should be made of the IV grade dental plaster, e.g. Student IV, which will ensure perfect reproduction and fitting in the patient's mouth.

Before starting the design of the denture, the paraleometer analysis should be performed, in order to determine the insertion path and examine the depth of the undercut on the model. Thanks to the flexibility of the material, it is possible to perform works using the undercuts of the gum tissue and the abutment tooth for retention. When designing the retention elements in the front section, the minimal tooth coverage is recommended for the sake of aesthetics. In the case of lateral teeth, the use of larger clasps will allow for greater stability of the denture. In order to achieve satisfactory retention, the undercuts of the gum tissue and teeth should be used. It is most convenient to use the pencil to apply the design to the model.

Preparation of the model for reproduction: The denture insertion path should be determined on the basis of the paraleometer analysis. It is necessary to block only those undercuts that cause difficulties when placing the denture on the prosthetic area. It is forbidden to block the undercuts on the tongue surface and the last teeth in the arch, because they can additionally improve the maintenance of the denture. It is necessary to relieve the necks of the teeth and sensitive areas with wax.

Duplication of the model: The model should be duplicates according to the standard procedures used when duplicating in the flask.

Teeth preparation: Use classic acrylic and composite teeth. Artificial teeth do not chemically bond to the thermoplastic material.

In order to properly maintain the teeth in the material, retention holes in the shape of the letter "T" should be drilled - 0.5 mm.

Artificial teeth do not chemically bond to the thermoplastic material. Corrections and repairs of the denture made of the Villacyl Thermo Press does not chemically bond with traditional acrylic materials.

Repairs and relining of polyamide-based thermoplastics. In order to perform the material injection in the case of repair or relining, apply a small amount of material for polyamide-based repairs before performing the injection to the previously developed edges of the tooth in the thermoplastic material.

Wax modeling: Before modeling, isolate the model using the UV isolator Villacyl Thermo Iso-Clear, and then cure in the UV polymerization lamp. The recommended minimum thickness of the modeled element is 1mm. The flexibility of the denture base elements can be adjusted by the thickness and shape of the wax modeling. The right design and optimal denture preparation process ensures that the denture made of the Villacyl Thermo Press will not break. The denture plate in the neck of the tooth should be modeled on the border of the acrylic teeth preparation. Do not model to precipitate the material.

Flasking: In the lower part of the flask, class I (Student III) or class IV (Student IV) stone should be used. To pour the countermold, it is necessary to use the class IV stone (Student IV). In order to facilitate the unflasking, apply the vaseline to the inner surface of the flask. Place the model in the central part of the flask, closer to the injection hole of the can. Care should be taken not to smudge the edge of the can. Wait for the channels to completely set. Apply U-shaped injection channels with a diameter of 5 mm to the denture. It is recommended that the channels adhere to the denture plate with the surface corresponding to the diameter. Apply the stone-stone isolate to the stone parts. Then assemble the flask, block and screw the flask halves using screws. Do not tighten the screws. Then fill the other part of the flask with the super hard stone and leave it to completely harden the stone.

Recommended location of injection channels (5mm):


Can scalding: The wax scalding procedure is the same as in the case of traditional acrylic dentures. Take special care when scalding the denture. Use protective clothing, gloves and goggles when working with hot water and hot can.

Heat the can in hot water or steam bath for 6-10 minutes. Pay attention not to exceed the scalding time, because the liquefied wax penetrates the plaster structure, which makes it difficult to isolate. After opening the flat, remove the wax residue thoroughly, using the stream of water from the scalding machine. Dry the flask from the water and cool the flask. Then, in order to obtain a smooth surface, it is necessary to apply the Villacyl Thermo Iso-Clear isolator to the countermold surface and injection channels, and then polymerize. Be careful not to insulate acrylic teeth. After closing and screwing the can with screws, heat the flask in the scalding machine to 95°C-100°C.

Material injection: Take special care when injecting the material. Use protective clothing, gloves and goggles when working on the injection molding machine. Failure to follow the recommendations may cause burns.

Injection parameters for the Villacyl Thermo Press: Observe the following parameters to achieve the optimal properties of the denture. However, they may differ depending on the efficiency and type of the injection machine used.

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Storage and handling:

The product is sensitive to UV radiation. Store the product in the original, airtight packaging and in dark and well-ventilated room at a

Table 1

Technical parameters for injection	
Injection temperature [°C]	275 - 285
Melting time [min]	15 - 27
Small cartridge (max. 13g)	15
Medium cartridge (max. 17g)	19
Big cartridge (max. 25g)	27
Injection pressure [bar]	7.5 - 9.0
Cooling time [min]	min. 15

The material injection procedure recommended by the manufacturer of the device in which the material will be injected, must be strictly observed!

Note! Place the can in the injection molding machine at a temperature from 90°C to 100°C.

Note! Do not set the injection temperature higher than 285°C.

Unflasking: After the process, let the flask cool down for about 30 minutes. Do not speed up the cooling process with water or other coolants.

Note! Do not hit the flask, not to damage its walls during unflasking. The damaged flask is not suitable for injection and cannot be used.

Before starting the design of the denture, the paraleometer analysis should be performed, in order to determine the insertion path and examine the depth of the undercut on the model. Thanks to the flexibility of the material, it is possible to perform works using the undercuts of the gum tissue and the abutment tooth for retention. When designing the retention elements in the front section, the minimal tooth coverage is recommended for the sake of aesthetics. In the case of lateral teeth, the use of larger clasps will allow for greater stability of the denture. In order to achieve satisfactory retention, the undercuts of the gum tissue and teeth should be used. It is most convenient to use the pencil to apply the design to the model.

Preparation of the model for reproduction: The denture insertion path should be determined on the basis of the paraleometer analysis. It is necessary to block only those undercuts that cause difficulties when placing the denture on the prosthetic area. It is forbidden to block the undercuts on the tongue surface and the last teeth in the arch, because they can additionally improve the maintenance of the denture. It is necessary to relieve the necks of the teeth and sensitive areas with wax.

Duplication of the model: The model should be duplicates according to the standard procedures used when duplicating in the flask.

Teeth preparation: Use classic acrylic and composite teeth. Artificial teeth do not chemically bond to the thermoplastic material.

In order to properly maintain the teeth in the material, retention holes in the shape of the letter "T" should be drilled - 0.5 mm.

Artificial teeth do not chemically bond to the thermoplastic material. Corrections and repairs of the denture made of the Villacyl Thermo Press does not chemically bond with traditional acrylic materials.

Repairs and relining of polyamide-based thermoplastics. In order to perform the material injection in the case of repair or relining, apply a small amount of material for polyamide-based repairs before performing the injection to the previously developed edges of the tooth in the thermoplastic material.

Wax modeling: Before modeling, isolate the model using the UV isolator Villacyl Thermo Iso-Clear, and then cure in the UV polymerization lamp. The recommended minimum thickness of the modeled element is 1mm. The flexibility of the denture base elements can be adjusted by the thickness and shape of the wax modeling. The right design and optimal denture preparation process ensures that the denture made of the Villacyl Thermo Press will not break. The denture plate in the neck of the tooth should be modeled on the border of the acrylic teeth preparation. Do not model to precipitate the material.

Flasking: In the lower part of the flask, class I (Student III) or class IV (Student IV) stone should be used. To pour the countermold, it is necessary to use the class IV stone (Student IV). In order to facilitate the unflasking, apply the vaseline to the inner surface of the flask. Place the model in the central part of the flask, closer to the injection hole of the can. Care should be taken not to smudge the edge of the can. Wait for the channels to completely set. Apply U-shaped injection channels with a diameter of 5 mm to the denture. It is recommended that the channels adhere to the denture plate with the surface corresponding to the diameter. Apply the stone-stone isolate to the stone parts. Then assemble the flask, block and screw the flask halves using screws. Do not tighten the screws. Then fill the other part of the flask with the super hard stone and leave it to completely harden the stone.

Recommended location of injection channels (5mm):

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