making work easy


- DE
- EN
- FR
- IT
- ES
- PT
- TR
- RU
- PL
- ZH
- JA
- KO
- AR
- UK

# SILENT TC

Made in Germany

21-6812 26112019

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making work easy

# SILENT TC

TRANSLATION OF THE ORIGINAL INSTRUCTIONS FOR USE

Made in Germany

21-6812 26112019

# Contents

EN

<b>1</b>	<b>Introduction</b>	<b>3</b>
1.1	Symbols used	3
<b>2</b>	<b>Safety</b>	<b>3</b>
2.1	Intended Use	3
2.2	Improper Use	3
2.3	Ambient Conditions for Safe Operation	4
2.4	Ambient Conditions for Storage and Transport	4
2.5	Hazard and Warning Information	4
2.5.1	General Information	4
2.5.2	Specific Information	5
2.6	Authorized Persons	5
2.7	Disclaimer	5
<b>3</b>	<b>Product Description</b>	<b>6</b>
3.1	General Description	6
3.2	Component Groups and Functional Elements	6
3.3	Delivery	7
3.4	Accessories	7
<b>4</b>	<b>Initial Operation</b>	<b>8</b>
4.1	Unpacking	8
4.2	Setting-up	8
4.3	Electrical connection	8
4.4	Connection to a Suction Point	8
4.5	Connecting the Electrical Devices	9
4.6	External Exhaust Air Duct	9
<b>5</b>	<b>Operation</b>	<b>9</b>
5.1	Switching On	9
5.2	Select operating mode: Automatic operation / continuous operation	9
5.3	Suction Power	10
5.4	Filter Cleaning	10
5.4.1	Automatic Filter Cleaning	10
5.4.2	Manual Filter Cleaning	10
5.5	Setting the Parameter	10
5.5.1	Acoustic signal (Buzzer)	11
5.5.2	Time Interval for "Full" Display Dust Drawer	11
5.5.3	Automatic Switch-on for Automatic Operation	12
5.5.4	After-run Time of the Suction Turbine	12
<b>6</b>	<b>Cleaning / Maintenance</b>	<b>13</b>
6.1	Cleaning	13
6.2	Empty Dust Drawer	13
6.3	Change Fine Filter	13
6.4	Self-diagnosis	14
6.5	Fuses	14
6.6	Spare Parts	14
6.7	Factory Settings	14
<b>7</b>	<b>Troubleshooting</b>	<b>15</b>
<b>8</b>	<b>Technical Data</b>	<b>16</b>
<b>9</b>	<b>Warranty</b>	<b>17</b>
<b>10</b>	<b>Disposal Information</b>	<b>17</b>
10.1	Disposing of Consumables	17
10.2	Disposing of the Unit	17
10.2.1	Disposal Instructions for countries in the EU	17

# 1 Introduction

## 1.1 Symbols used

In the instructions for use and on the unit itself you will find these symbols with the following meanings:

### **Danger**

**Immediate risk of injury. Consult accompanying documents!**

### **Electrical current**

**This indicates a risk of injury due to electrical current.**

### **Caution**

**Risk of damage to the unit if the instruction is not followed.**

### **Notice**

**A general notice that improves and eases use.**

**The appliance complies with applicable EU directives.**

**Within the EU, this device is subject to the provisions of the directive 2002/96/EC (WEEE directive).**

### ► **List, particular attention should be paid**

- List
- List

⇒ Instructions / appropriate action / input / operational sequence:

You will be asked to carry out the action in a specified order.

- ◆ Result of an action / reaction of the device / reaction of the program:

The unit or program reacts as a result of your actions or when a specific incident occurs.

Other symbols are explained as they occur.

## 2 Safety

### 2.1 Intended Use

This device is designed to extract dry, non-explosive dust.

The unit is intended solely for use in a commercial dental laboratory and dental practice.

The intended use also includes compliance with the instructions specified by the manufacturer concerning operation, servicing and maintenance.

### 2.2 Improper Use

**The equipment is not suitable for use with dental CAM units!**

**For information concerning SILENT extraction units for CAM machines, please refer to [www.renfert.com](http://www.renfert.com).**

Fire-promoting, easily flammable, red-hot, burning or explosive materials must not be suctioned into the device.

It is not permitted to suction liquids.

This device is not intended for private, household use.

Any use other than specified in these instructions is deemed improper and constitutes a misuse of the device.

The manufacturer shall not be liable for damages caused by improper use.

Only spare parts and accessories supplied or authorized by Renfert GmbH may be used with this product. If other spare parts or accessories are used, this could have a detrimental effect on the safety of the device, increase the risk of serious injury and lead to damage to the environment or the device itself.

## 2.3 Ambient Conditions for Safe Operation

The device may only be operated:

- Indoors,
- Up to an altitude of 2,000 m above sea level,
- At an ambient temperature of between 5 - 40 °C [41 - 104 °F] \*),
- at a maximum relative humidity of 80 % at 31 °C [87.8 °F], linear reduction up to 50 % relative humidity at 40 °C [104 °F] \*),
- With mains power where the voltage fluctuations do not exceed 10 % of the nominal value,
- Under contamination level 2 conditions,
- Under over-voltage category II conditions.

\*) Between 5 - 30 °C [41 - 86 °F] the device can be operated at a relative humidity of up to 80 %. At temperatures between 31 - 40 °C [87.8 - 104 °F] the humidity must decrease proportionally in order to ensure operational readiness (e.g. at 35 °C [95 °F] = 65 % humidity, at 40 °C [104 °F] = 50 % humidity). The device may not be operated at temperatures above 40 °C [104 °F].

## 2.4 Ambient Conditions for Storage and Transport

For storage and transport the following specifications to ambient conditions apply:

- Ambient temperature - 20 – + 60 °C [- 4 – + 140 °F].
- Maximum relative humidity 80 %.

## 2.5 Hazard and Warning Information

### 2.5.1 General Information

- ▶ If the device is not used in compliance with the supplied instructions, the safety of the device can no longer be guaranteed.
- ▶ The device may only be operated using a mains cable with the country-specific plug system. Any necessary alterations must be carried out by a qualified electrician.
- ▶ The device may only be operated if the information on the identification plate conforms to the specifications of your local mains power supply.  
After removing the Dust Drawer, you can find the identification plate below, inside the device on the left side.
- ▶ The device may only be plugged into outlets which are connected to the protective conductor system.
- ▶ The mains plug must be easily accessible.
- ▶ Disconnect the device from the mains before carrying out work on the electrical parts.
- ▶ Check connection cables (such as power supply cords), tubes and housing (i.e. the key-pad) regularly for damage (i.e. kinks, cracks and porosity) or signs of ageing. Do not operate units with damaged connection cables, hoses or housing parts or other defects!
- ▶ Defective devices must be put out of service immediately. Remove the mains plug and ensure the device is not used. Send the device in for repair!
- ▶ Only operate the device under supervision.
- ▶ Please observe the national accident prevention regulations!
- ▶ It is the responsibility of the operator that national regulations during operation and regarding a repeated safety inspection of electrical equipment are complied with. For Germany these are the regulation 3 by DGUV (German Statutory Accident Insurance) in relation with VDE 0701-0702 (Association for Electrical, Electronic and Information Technology).
- ▶ Information on REACH and SVHC is available on our website [www.renfert.com](http://www.renfert.com), in the Support area.

## 2.5.2 Specific Information

- ▶ The mains socket on the device is only designed for purposes specified in the operating instructions. Connecting other devices may cause material damage.
- ▶ Switch off the electrical device before connecting an electrical appliance to the appliance outlet.
- ▶ Read the operating instructions for the electrical appliance which is to be connected and observe the stipulated safety instructions.
- ▶ Please observe the national regulations and permitted exposure to dust in a working environment. Please ask the “National Institute for Occupational Safety and Health” or other responsible authority.
- ▶ Always refer to the relevant safety data sheets, when extracting hazardous materials.
- ▶ Always wear protective gear, when extracting hazardous materials.
- ▶ It is necessary to wear suitable personal protective equipment when emptying the dust drawer or cleaning, depending on the type of extracted material.
- ▶ When disposing of the extracted material or used filter, please observe the local specifications and accident prevention regulations!
- ▶ Make sure the dust drawer is fully closed during operation.
- ▶ Do not operate without a suction hose.
- ▶ Do not extract flammable or explosive gasses, fumes or dust.
- ▶ Do not extract hot materials.
- ▶ Do not extract liquids.
- ▶ If the dust extractor is employed to extract hazardous materials, appropriate personal protective gear must be worn and steps must be taken to ensure that the exhaust air is properly ventilated. Please refer to the associated safety data sheets for specific requirements.
- ▶ Dispose of extracted material according to local statutory regulations.

## 2.6 Authorized Persons

Operation and maintenance of the device may only be performed by qualified personnel.

Minors and pregnant women may only operate and service the device if they are wearing appropriate protective gear, in particular if the device is being used to extract hazardous materials.

Any repairs not specifically described in these operating instructions may only be carried out by a qualified electrician.

## 2.7 Disclaimer

Renfert GmbH shall be absolved from all claims for damages or warranty if:

- ▶ The product is employed for any purposes other than those specified in the operating instructions.
- ▶ The product is altered in any way other than those alterations described in the operating instructions.
- ▶ The product is not repaired by an authorized facility or if non-original Renfert parts are implemented.
- ▶ The product continues to be used despite obvious safety faults or damage.
- ▶ The product is subjected to mechanical impacts or is dropped.

# 3 Product Description

## 3.1 General Description

This is a workplace and appliance suction unit for extracting dust in dental laboratories. The extraction unit can be operated both manually and automatically, depending on the connected, electrical, dust-generating equipment.

## 3.2 Component Groups and Functional Elements

- 1 SILENT TC
- 2 Control panel
- 3 On / off switch
- 4 Dust drawer
- 4a Clamp fastener
- 5 Fine filter
- 6a Suction port
- 7 Velcro® strip
- 8 Device socket
- 9 Mains cable
- 10 Protective device switch (2 x)
- 11 Exhaust air filter / Exhaust air outlet
- 12 Suction hose


Fig. 1


Fig. 2

- |   |  |
|---|--|
| <b>3</b> On / off switch  | <b>24</b> Display                                  |
| <b>20</b> Display Automatic operation                           | <b>25</b> [+] -key (increase value)                |
| <b>21</b> Start operation key, Automatic / continuous operation | <b>26</b> Display empty dust drawer                |
| <b>22</b> Display continuous operation                          | <b>27</b> Enter key, save input                    |
| <b>23</b> [-] key (reduce value)                                | <b>28</b> Cover (do not remove, service interface) |

### 3.3 Delivery

- 1 SILENT TC
- 1 Quick Start Guide
- 1 Suction hose, 3 m, incl. 2 end mufflers

### 3.4 Accessories

- 2934 0014 Waste bag for TC, TC2, EC, EC2 & PowerCAM EC (25 pcs.)
- 2934 0015 Waste bag for SILENT TC, TC2, EC & EC2 Starter Kit
- 2921 0003 End bushing set, 2 pieces
- 90003 4240 Suction hose, 3 m, incl. 2 end mufflers
- 90003 4826 Suction hose, antistatic, 3 m, incl. 2 end mufflers
- 90115 0823 Suction hose inside diameter 38 mm, 6 m
- 90215 0823 Suction hose inside diameter 38 mm, 9 m
- 90003 4305 Adapter for hose connection
- 90003 4430 Universal suction hose adapter
- 90003 4314 Y-junction
- 2925 0000 Extractor clamp
- 2925 1000 Glass pane with attachment (for the suction hood)
- 2926 0000 Y-Junction for dust extractor
- 2937 0002 External exhaust air vent for SILENT TC/EC/PowerCAM
- 2934 0007 90° angled connector Silent

For more details or further accessories, please see [www.renfert.com](http://www.renfert.com).

## 4 Initial Operation

### 4.1 Unpacking

- ⇒ Remove the device and accessories from the package.
- ⇒ Check the delivery for completeness (compare the list included in delivery).

### 4.2 Setting-up

The suction unit is a standing device which must not be operated in a horizontal position.

Place the suction unit so that:

- the exhaust air outlet (11, Fig. 1) is not blocked.
- the front side is easily accessible for removing the dust drawer.

**If the suction unit is placed in a closed cabinet, the warm exhaust air must be transported out of the cabinet using one of the measures shown in Figs. A, B, C at the start of this manual.**

- External exhaust air vent (A) (see chapter 4.6).
- Opening in the cabinet wall (B), min. 250 x 120 mm, directly opposite the exhaust air outlet (11, Fig. 1).
  - Distance of the cabinet back side to the wall: min. 100 mm,
  - Distance from the extraction unit to the cupboard back side: max. 25 mm.
- Remove the back side of the cabinet (C), distance from the back of the cabinet to the wall min. 50 mm.

If the warm exhaust air is transported out through the openings to the rear, it is essential that the air can escape unhindered from there.

### 4.3 Electrical connection

**Before making the electrical connection, check that the voltage indicated on the rating plate matches the local voltage supply.**

**Ensure that the layout of live parts (socket, plug and coupling) and extension cables are arranged so that the protection class is maintained.**

- ⇒ Switch the device off at the on / off switch (3, Fig. 2).
- ⇒ Unroll the mains cable (9, Fig. 1) and insert the plug into the building's electrical socket.

### 4.4 Connection to a Suction Point

**Warning risk of injury!**  
**When shortening the suction tube, please be sure to cut the integrated wire as straight as possible.**

**Long suction tubes, tight bends and kinks greatly reduce the suction power at the suction point.**

- ⇒ Shorten the suction tube if possible.
- ⇒ Open the adhesive strap (7).
- ⇒ Plug the suction tube (12) into the suction socket (6).
- ⇒ Fix the suction tube in place with the adhesive strap.
- ⇒ Connect the suction tube to the correct suction point.

**If the diameter is not suitable, use an adapter (see accessories), to avoid loss of suction.**

**Void steep sloping and "sagging" of the tubes.**


Fig. 3

## 4.5 Connecting the Electrical Devices

⇒ Connect the electronic device to the appliance socket (8) on the back.

**When connecting electrical appliances to the suction unit, ensure that the maximum permissible power for connected devices is not exceeded (see chapter 8, Technical data).**


Fig. 4

## 4.6 External Exhaust Air Duct

The exhaust air can be transported out of the laboratory via an external exhaust air duct (see accessories).

The installation instructions are included in the external exhaust air duct.

**When using suction devices in conjunction with an external exhaust air duct, a large amount of air is extracted from the room per hour.**

**This can lead to a negative pressure, whereby, when using air-dependent fireplaces with gas, liquid or solid fuels, poisonous gases (e.g. carbon monoxide) will be sucked into the (working) space.**

**Depending on the structural situation, it is therefore essential to ensure there is an additional air supply or negative pressure monitoring which is checked by the appropriate institutions (e.g. chimney sweep or heating engineer).**

# 5 Operation

The suction is controlled by the keys on the control panel (Fig. 2).

## 5.1 Switching On

The suction unit is switched on and off via the on / off switch (3).

- ◆ The suction unit performs an automatic filter cleaning sequence (for approx. 8 seconds, loud vibrating noise). (see chapter 5.4).

Afterwards, the suction unit is in the last set mode.


Fig. 5

## 5.2 Select operating mode: Automatic operation / continuous operation

The extraction unit has two operating modes.

The set operating mode is shown in the display (20) / (22).

- Automatic mode (20):  
The suction is dependent on the electrical device connected to the appliance outlet (8, Fig. 1).
- Permanent operation (22):  
Suction running continuously.

⇒ Press operating mode key (21).

- ◆ Changing the operating mode.


Fig. 6

**To start or stop the unit during permanent operation, use the operating mode key (21). Do not use the on / off switch for this.**

**If the suction to a connected electrical appliance (e.g. a handpiece) does not react as required (e.g. does not activate, even though a connected device is in operation) then the switch-on threshold for automatic operation must be adjusted (see chapter 5.5.3).**

## 5.3 Suction Power

The suction power for the suction unit can be set in 4 levels.

The current level is shown in the display (24).

Change the suction power:

⇒ [-] Press key (23)

- ◆ Reduce suction

⇒ [+] Press key (25)

- ◆ Increase suction

The suction always starts with the last used suction level.


Fig. 7

## 5.4 Filter Cleaning

In order to ensure maximum suction performance, the suction unit has a device for cleaning the fine filter unit. The cleaning takes approx. 8 sec.

The cleaning is carried out:

- automatic
  - After the unit is switched on;
  - With insufficient suction power (flow velocity falls below an internal limit value);
  - If the unit has been operated for more than 8 hours (turbine running time) without switching it off;
- Manually, e.g. before removing the dust drawer to empty it (see chapter 6.2).

### 5.4.1 Automatic Filter Cleaning

⇒ Limit value is too low:

- ◆ Suction turbine is stopped.
- ◆ Display dust drawer (26) blinks.
- ◆ The outer segments of the display (24) light up in a circulating manner
- ◆ An acoustic signal indicates the filter cleaning is about to commence.
- ◆ The filter cleaning function is performed.
- ◆ Suction turbine starts again.

### 5.4.2 Manual Filter Cleaning

⇒ Press Enter key (27) for 2 seconds.

- ◆ The filter cleaning function is performed.

## 5.5 Setting the Parameter

On delivery, there are parameters set, which in most cases, make it possible to work with the suction unit without problems. Only in individual cases, if this is not possible, the parameters can be changed.

The setting of various parameters and the execution of a self-diagnosis are carried out in the programming mode.

The display (24, Fig. 2) shows which parameter are set:

	"d"	Perform self-diagnosis (see chapter 6.4).
	"b"	Acoustic signal on / off (decimal point on / off).
	"A"	Time Interval for "Full" Display Dust Drawer.
	"c"	Setting the automatic switch-on (calibration).
	"t"	Set after-run (time) for the suction turbine.

In order to set the different parameters, the programming mode must be started and the parameter selected.

Confirm and save settings with the Enter key (27). Successful saving is acknowledged with an acoustic signal.

If a change is not to be applied, press the mode key (21) again to cancel.


Fig. 8

### 5.5.1 Acoustic signal (Buzzer)

When on, various inputs are confirmed with an acoustic signal.

To change:

- ⇒ Press operating mode key (21) for 3 sec.
  - ◆ Starting the programming mode.
- ⇒ Press operating mode key 1 x again.
  - ◆ In the display a „b“ blinks.
- ⇒ Press Enter key (27).
  - ◆ Acoustic Signal selected.
  - ◆ The decimal point blinks: The acoustic signal is switched on.
  - ◆ The decimal point is off: The acoustic signal is switched off.
- ⇒ Switch the acoustic signal on or off with the [+] / [-] keys.
- ⇒ Press Enter key (27).
  - ◆ Save changes, exit programming mode.


Fig. 9

### 5.5.2 Time Interval for “Full” Display Dust Drawer

After reaching a set time interval, a request appears to empty the dust drawer.

It is possible to select one of 5 different times.

Time interval / hrs	Value in display (24, Fig. 2)
2	1
5	2
10	3
50	4
100	5

Factory setting

***In the case of heavy dust generators (e.g. blasting equipment), the time interval “emptying the dust drawer” should be set to 5 hours, if necessary to 10 hours.***

- ⇒ Press operating mode key (21) for 3 sec.
  - ◆ Starting the programming mode
- ⇒ Press operating mode key 2 x again.
  - ◆ In the display “A” blinks.
- ⇒ Press Enter key (27).
  - ◆ Time Interval for “Full” display Dust Drawer selected.
  - ◆ In the display (24), “A” and a number for the time interval will be shown alternately.
- ⇒ Set the required time interval with the [+] / [-] keys.
- ⇒ Press Enter key (27).
  - ◆ Save changes, exit programming mode.


Fig. 10

### 5.5.3 Automatic Switch-on for Automatic Operation

Only delivery, the switch-on automatic is set in such a way that operation with the most handpieces is possible.

If the suction to a connected electrical appliance (e.g. a handpiece) does not react as required (e.g. does not activate, even though a connected device is in operation) then the switch-on threshold for automatic operation must be adjusted.

⇒ Press operating mode key (21) for 3 sec.

- ◆ Starting the programming mode

⇒ Press operating mode key 3 x again.

- ◆ In the display “c” blinks

⇒ Press Enter key (27).

- ◆ Setting the automatic switch-on.

- ◆ The display will flash alternately “c” (calibration) and “0”.

⇒ Switch off electrical devices without standby mode.

⇒ For electrical devices with a standby mode (e.g. handpieces), switch to standby mode (e.g. with a handpiece, switch only the control unit on without activating the handpiece).

⇒ Press the Enter key (27).

- ◆ The display will flash alternately “c” (calibration) and “l”.

⇒ An electrical device without standby mode must be switched on.

⇒ Operate the handpiece at the speed at which the suction is to be started and run for approx. 3 - 5 seconds (so that overcurrent during start-up is not detected).

⇒ Press the Enter key (27) (whilst the handpiece is still running).

- ◆ An acoustic signal acknowledges the information.
- ◆ The setting is saved.
- ◆ Exit programming mode.

**For devices with standby mode, these should be switched on approx. 5 minutes before setting to avoid incorrect measurements.**

**If there is an error with the switch-on function, the display “C” will blink and an acoustic signal is emitted twice. Set the changes again.**

### 5.5.4 After-run Time of the Suction Turbine

After-run time = the time between switching off a connected consumer and turning off the suction turbine.

You can select between 5 times (factory setting 3 sec.).

After-run time / sec.	Value in display (24, Fig. 2)
0	1
3	2
5	3
10	4
20	5
30	6

Factory setting

To change:

⇒ Press operating mode key (21) for 3 sec.

- ◆ Starting the programming mode.

⇒ Press the operating mode key 4 x.

- ◆ In the display “t” blinks.

⇒ Press Enter key (27).

- ◆ Setting the after-run time is selected.

- ◆ In the display (24), “t” and a number for the time interval will be shown alternately.

⇒ Set the required time interval with the [+] / [-] keys.

⇒ Press Enter key (27).

- ◆ Save changes, exit programming mode.


Fig. 12

**In order to prevent the suction from being switched on / off too often in the case of short interruptions (for example when working with the handpiece), a higher time value can be selected.**

## 6 Cleaning / Maintenance

The device does not contain any internal parts that require maintenance. Opening the device other than described in the instructions is not permitted!

### 6.1 Cleaning

To clean the device wipe only with a damp cloth.  
Do not use and solvent-based or abrasive cleaning agents.

### 6.2 Empty Dust Drawer

After reaching the set time interval (see Chapter 5.5.2), the dust drawer is requested to be emptied:

- ◆ An acoustic signal sounds 3 x.
- ◆ Display dust drawer (26, Fig. 2) lights up.
- ◆ The display (24, Fig. 2) shows a container filling up.

Before removing the dust drawer, a filter cleaning function should be performed:

- ⇒ Press Enter key (27, Fig. 2) for 2 seconds.
  - ◆ The filter cleaning function is performed.
  - ◆ The counter for recording the time interval is reset.

After the filter cleaning has finished:

- ⇒ Unlock the dust drawer with the clamp fastener (4a).
- ⇒ Pull the dust drawer (4) forward and empty.
- ⇒ Replace the dust drawer again until it locks into place.
- ⇒ Close the clamp fastener (4a).

If the dust drawer has been emptied without prior filter cleaning, the counter must be reset separately to record the time interval.

- ⇒ Press the Enter key (27, Fig. 2) (counter for recording the time interval is reset).
  - ◆ Acoustic signal acknowledges the information.
  - ◆ Dust drawer display (26, Fig. 2) is deleted.


Fig. 13

**If the dust drawer is not emptied, the dust drawer (26, Fig. 2) display will continue to light up. After switching the device off and on again, a 3 x acoustic signal is used to indicate the emptying of the dust drawer.**

### 6.3 Change Fine Filter

The effectiveness of the filter cleaning is determined by monitoring the flow rate (internal pressure sensor).

If a filter cleaning is required several times in succession after less than 2 hours, this means that the fine filter is so highly clogged that the filter cleaning does not show sufficient effect and the fine filter has to be changed.

This is indicated as follows:

- ◆ "F" in display (24, Fig. 2).
- ◆ For a period of 15 min., 2 long beeps are sounded for 3 min at a time.

The error message can be switched off by switching off the suction.

**If requested to change the fine filter or the suction performance does not improve substantially despite repeated cleaning or only for a short period, the fine filter must be changed soon. Further operation can lead to damage to the device.**

When replacing the fine filter, pay attention to the correct fit, otherwise leaks will occur. Please refer to the assembly instructions at the end of the manual, which are also included with the new fine filter.

## 6.4 Self-diagnosis

With the help of the self-diagnosis, the control system checks the function of the turbine and parts of the electronics.

If no fault is detected, the suction is ready for operation after the self-diagnosis has finished. If an error is detected, this is indicated.

Self-diagnosis can be started manually and then performed once.

⇒ Press operating mode key (21) for 3 sec.

- ◆ In the display "d" blinks.

⇒ Press Enter key (27).

During self-diagnosis:

- An acoustic signal sounds.
- All displays are switched on shortly.
- A "d" (Diagnosis) is shown in the display (24, Fig. 2).
- The control program checks various internal components.
- The suction turbine is switched on shortly.


Fig. 14

## 6.5 Fuses

The protection of the suction device is achieved by two device protection switches (10, Fig. 1).

A triggered device protection switch is reset by depressing the button.

**Repeated triggering of a device protection switch is caused by a defect within the device. Send the device in for repair!**

## 6.6 Spare Parts

You can find components subject to wear and the spare parts on the spare part list in the internet at [www.renfert.com/p918](http://www.renfert.com/p918)

Enter the following item number: 29350000.

The components excluded from the warranty (such as consumables or parts subject to wear and tear) are marked on the spare part list.

Serial number and date of manufacturing are shown on the type plate of the unit.

## 6.7 Factory Settings

⇒ Switch the appliance off.

⇒ Press the [-] key and [+] key at the same time.

⇒ Switch device on and press the keys for 3 seconds.

- ◆ All 4 displays blink 3 x.
- ◆ All values are reset to the factory settings.
- ◆ The suction unit performs an automatic filter cleaning.

**Factory settings:**

Function / Characteristic	Setting range	Factory setting
Operating mode	Automatic mode / Permanent operation	Automatic mode
Suction level	1 - 4	2
Time Interval for Dust Drawer	2 - 100 hrs	50 hrs
After-run time	0 - 30 sec.	3 sec.
Automatic switch-on for automatic operation	1 - 100 W	8 W

# 7 Troubleshooting

Trouble	Reason	Solution
<b>When the device is switched on, there is a loud vibration noise for approx. 8 seconds.</b>	<ul style="list-style-type: none"> <li>The suction unit performs an automatic filter cleaning sequence.</li> </ul>	<ul style="list-style-type: none"> <li>This is functional and requires no intervention.</li> </ul>
<b>An acoustic signal is emitted, the suction unit is switched off and a filter cleaning is performed.</b>	<ul style="list-style-type: none"> <li>The internal flow velocity limit was not reached.</li> <li>After 8 hours of operation (turbine air time) without switching off, a filter cleaning is carried out.</li> </ul>	<ul style="list-style-type: none"> <li>Repeat work after cleaning is finished.</li> <li>Switch off devices at the end of each working day at the on / off switch (3, Fig. 1).</li> </ul>
<b>After switching on, the display shows "Dust drawer" (26, Fig. 2) and a 3 x acoustic signal is heard.</b>	<ul style="list-style-type: none"> <li>The time interval for emptying the dust drawer has expired and the dust drawer has not yet been emptied.</li> <li>The emptying of the dust drawer has not yet been confirmed.</li> </ul>	<ul style="list-style-type: none"> <li>Empty Dust Drawer and confirm by pressing the Enter key (27, Fig. 2).</li> <li>Confirm emptying the Dust Drawer by pressing the Enter key (27, Fig. 2).</li> </ul>
<b>The display shows "C".</b>	<ul style="list-style-type: none"> <li>The automatic switch-on function cannot be set. At the selected rotational speed on the handpiece, the difference between stand-by current and operating current is too low.</li> <li>The automatic switch-on function cannot be set.</li> </ul>	<ul style="list-style-type: none"> <li>Adjust the switch-on threshold to a higher speed on the handpiece.</li> <li>Device in permanent operation.</li> </ul>
<b>The display shows "H".</b>	<ul style="list-style-type: none"> <li>The electronics have become overheated.</li> </ul>	<ul style="list-style-type: none"> <li>Switch off the device and let it cool down.</li> <li>Provide adequate cooling, e.g.                             <ul style="list-style-type: none"> <li>Observe device set-up in chapter 4.2.</li> <li>Use the external exhaust function (see chapter 4.6).</li> <li>Change the Fine Filter (see Accessories and chapter 6.3).</li> </ul> </li> </ul>
<b>The display shows "E".</b>	<ul style="list-style-type: none"> <li>Leakage: Dust drawer not airtight.</li> <li>Faults in electronics.</li> <li>The suction turbine is not running.</li> </ul>	<ul style="list-style-type: none"> <li>Check that the Dust Drawer fits correctly (see chapter 6.2).</li> <li>Contact Renfert / Service.</li> <li>When the suction turbine has been changed, check the plug-in contact at the suction turbine.</li> </ul>
<b>The display shows "F".</b>	<ul style="list-style-type: none"> <li>The fine filter is so heavily clogged that the filter cleaning does not show sufficient effect.</li> </ul>	<ul style="list-style-type: none"> <li>Change Fine Filter (see chapter 6.3).</li> </ul>
<b>The suction performance is insufficient.</b>	<ul style="list-style-type: none"> <li>Suction performance is set too low.</li> <li>Blockage or leakage in the Suction tube.</li> <li>Dust Drawer not airtight.</li> <li>Fine filter blocked.</li> </ul>	<ul style="list-style-type: none"> <li>Select a higher suction level.</li> <li>Check suction tube.</li> <li>Please note the information in chapter 4.4.</li> <li>Check the correct fit of the Dust Drawer (see chapter 6.2).</li> <li>Clean the filter.</li> <li>Switch the unit off and then on again, so that the filter cleaning process can be carried out. Change the fine filter (see chap. 6.3) (if the filter cleaning function does not improve the suction performance).</li> </ul>
<b>Dust Drawer over full.</b>	<ul style="list-style-type: none"> <li>Time interval to the "Empty Dust Drawer" has been set too high.</li> </ul>	<ul style="list-style-type: none"> <li>Set a shorter time interval (see chapter 5.5.2).</li> </ul>
<b>The signal for emptying the dust drawer appears, even though it is not yet full.</b>	<ul style="list-style-type: none"> <li>The set time interval to the "Empty Dust Drawer" is too low.</li> </ul>	<ul style="list-style-type: none"> <li>Set a larger time interval (see chapter 5.5.2).</li> </ul>

Trouble	Reason	Solution
Suction starts even though the electrical appliance connected to the socket is not being used.	• The automatic switch-on function has been set too low.	• Setting the automatic switch-on (see chapter 5.5.3).
Suction does not stop when the electrical appliance is switched off.	• The automatic switch-on function has been set too low.	• Setting the automatic switch-on (see chapter 5.5.3).
Suction does not start even though the electrical appliance connected to the socket is in used.	• The automatic switch-on function has been set too high.	• Setting the automatic switch-on (see chapter 5.5.3).
Fuse triggers when a connected electrical device is switched on.	• Power consumption of the connected device too high.	• Observe the maximum connection value (see chapter 8).
The suction stops abruptly in automatic or continuous operation and the associated displays (20 / 22, Fig. 2) are still on.	• Suction turbine is overheated.  • Suction turbine is defect.	• Switch off the device and leave to cool for at least 60 min. • Check whether suction tube is blocked, eliminate blockage. • Switch the device off and then on again to perform a filter cleaning. Change fine filter (see chapter 6.3) if filter cleaning does not improve the performance.  • Change Suction turbine.
After the suction motor was changed, the suction unit stopped functioning or stops suddenly.	• The plug from the suction motor is not sufficiently clicked in place.	• Connect the suction motor plug correctly. Ensure that the plug is fitted correctly and clicked into place by pulling it.

## 8 Technical Data

Item no.	<b>2935 0000</b>	
Nominal voltage	<b>230 V</b>	
Permitted voltage:	220 - 240 V	
Mains frequency:	50 / 60 Hz	
Suction turbine power: *)	850 W	
Max. connection value of the appliance socket: *)	2000 W	
Total connected power: *)	2850 W	
Mains input fuse:	2 x 15 A(T)	
LpA * *) (at max. flow):	54,3 dB(A)	
Ø suction support: - inner - outer	45 mm [1.77 inch] 50 mm [1.97 inch]	
max. Volume flow:	3980 l/min [2.34 ft <sup>3</sup> /s]	
Vacuum pressure, max.:	281 hPa [4.1 psi]	
Fine filter: - Filter surface approx. - Filter quality	0.9 m <sup>2</sup> [1390 sq inch] Class M according to EN 60335-2-69	
Fill level dust drawer, approx.:	7 l [1.85 US gal]	
Weight (empty), approx.:	22,0 kg [48.5 lbs]	
Dimensions (width x height x depth):	275 x 535 x 540 mm [10.8 x 21.1 x 21.3 inch]	

\*) Performance values at nominal voltage

\*\*) Noise level according to EN ISO 11202

## 9 Warranty

Provided the unit is properly used, Renfert **warrants all components for 3 years.**

**The suction motor is guaranteed for 3 years, or a maximum of 1000 operating hours (engine running time).**

Warranty claims may only be made upon presentation of the original sales receipt from the authorized dealer. Parts which are subject to natural wear and tear (wear parts) and consumables are excluded from the guarantee. These parts are marked in the spare part list.

The warranty is voided in the case of improper use; failure to observe the operating, cleaning, maintenance and connection instructions; in case of independent repairs or repairs by unauthorized personnel; if spare parts from other manufacturers are employed, or in case of unusual influences or influences not in compliance with the utilization instructions.

Warranty service shall not extend the original warranty.

## 10 Disposal Information

### 10.1 Disposing of Consumables

Full dust bags and filters must be disposed of under compliance with locally applicable regulations. Depending on the material trapped by the filters, protective gear may need to be worn during disposal.

### 10.2 Disposing of the Unit

The unit must be disposed of by an authorized recycling operation. The selected firm must be informed of all possible health hazardous residues in the unit.

#### 10.2.1 Disposal Instructions for countries in the EU

To conserve and protect the environment, prevent environmental pollution and improve the recycling of raw materials, the European Commission adopted a directive that requires the manufacturer to accept the return of electrical and electronic units for proper disposal or recycling.

**Within the European Union, units with this symbol should therefore not be disposed of in unsorted domestic waste.**

Please contact your local authorities for more information on proper disposal.

**Hochaktuell und ausführlich auf ...  
Up to date and in detail at ...  
Actualisé et détaillé sous ...  
Aggiornato e dettagliato su ...  
La máxima actualidad y detalle en ...  
Актуально и подробно на ...**

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Renfert GmbH • Untere Gießwiesen 2 • 78247 Hilzingen/Germany